
5027 Sampler User Manual

5000 3694 / Rev. 3

5400029a

Information about revisions to this manual is sent out by Technical Notes.

Revision	Date of Issue	Revised Material
3	2005-04-21	Design modifications

NOTE: Please note that all information is liable to change without prior notice.

Please inform Technical Information Group at the address given below, if you have any opinions about or proposals for changes to this manual.

Technical Information Group

e-mail tig@foss.dk

FOSS Analytical AB, Box 70, SE-263 21 Höganäs, Sweden

Tel Int +46 42 361500, Fax +46 42 340349

English**Warning**

In order to find out the nature of the potential hazard, please consult this manual in all cases where this symbol is used.

The information will be found on the yellow pages.

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1	General	1:1
1.1	Introduction	1:1
1.2	Warranty Statement	1:1
1.3	Safety Precautions	1:1
2	Technical Description	2:1
2.1	General	2:1
2.2	Features	2:1
2.3	Technical Specification	2:1
3	Unpacking and Checking	3:1
4	Installation	4:1
4.1	General	4:1
4.2	Voltage Feed	4:1
4.3	Connecting the Sample Probe	4:2
4.4	Fitting the Sample Racks	4:2
4.5	Environmental Conditions	4:2
5	Operation	5:1
5.1	Connecting Analysis Systems	5:1
5.2	Identification of Operating Controls	5:1
5.3	Setting the Analysis Configuration	5:1
5.4	Starting the Sampler	5:2
5.5	Operation	5:2
6	Maintenance	6:1
6.1	Instructions for Cleaning	6:1
7	Accessories and Consumables	7:1
8	SAFETY INFORMATION (GB)	8:1
9	SICHERHEITSINFORMATIONEN (GE)	9:1
10	SIKKERHEDSINFORMATION (DK)	10:1
11	INFORMACION DE SEGURIDAD (ES)	11:1
12	TURVALLISUUSTIEDOT (FI)	12:1

13	CONSIGNES DE SECURITE (FR)	13:1
14	ΠΛΗΡΟΦΟΡΙΕΣ ΑΣΦΑΛΕΙΑΣ (GR)	14:1
15	UPPLÝSINGAR UM ÖRYGGI (IS)	15:1
16	INFORMAZIONI SULLA SICUREZZA (IT)	16:1
17	VEILIGHEIDSVOORSCHRIFTEN (NL)	17:1
18	SIKKERHETSINFORMASJON (NO)	18:1
19	INFORMAÇÕES DE SEGURANÇA (PT)	19:1
20	SÄKERHETSINFORMATION (SE)	20:1

1 General

1.1 Introduction

This manual contains material for the 5027 Sampler Unit.

1.2 Warranty Statement

In the agreements below, FOSS Analytical AB is mentioned as “the Seller”. Information in this manual is subject to change without prior notice.

Terms of Limited Warranty

All FOSS Analytical products are warranted against defects in materials and workmanship under normal use in accordance with FOSS Analytical terms of sale. The foregoing warrant shall not apply to defects resulting from improper or inadequate maintenance, adjustment, calibration or operation by Buyer.

The warranty period is one year for replacements parts. The warranty period becomes effective at installation or 3 months from date of shipment, whichever occurs first. No other warranty is expressed or implied.

FOSS Analytical AB further guarantees that the goods have been thoroughly checked and tested before shipment and that they are free of mechanical and electrical defects. Should defects due to faulty material or workmanship develop during the WARRANTY PERIOD, parts will be supplied on a no-charge basis in accordance with FOSS Analytical terms of sale (except for expendable materials).

In no event will FOSS Analytical AB or its agents be liable for any consequential, incidental or indirect damages arising out of the use or inability to use the 5027 Sampler Unit. The 5027 Sampler Unit should only be used as described in this manual.

Modifications, alterations, rebuilding or use of parts not authorized by FOSS Analytical AB violate the warranty. FOSS Analytical AB cannot be held responsible for damage to property or injury to persons occurring as a result of such use or actions.

The Seller accepts no contingent liability and shall in no event be liable for consequential damages.

Warranty Disclaimer: Local market regulations and market conditions may require a change of the above stated warranty. Please check with the local FOSS Analytical distributor to confirm the prevailing warranty policy for your area.

1.3 Safety Precautions

Warning

This device is equipped with a grounding/earthing type power plug for your protection against electrical shock hazard and should only be attached to a properly grounded/earthed receptacle.

Warning

Before replacing the fuses, disconnect incoming mains supply.

Warning

Electrical shock hazard. Covers or panel(s) should be removed by qualified personnel only.

Warning

The responsible body shall be made aware that, if the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Modifications, alterations, rebuilding or use of safety parts not authorized by FOSS Analytical AB violates the warranty. FOSS Analytical AB has no responsibility for damages, material or personal, occurring as a result of such actions.

2 Technical Description

2.1 General

The 5027 Sampler can be used for several of FOSS Analytical's analysis systems including SCABA, Aquatec®, FIAstar®, GrapeScan, MilkoScan and WineScan. The Sampler contains a microprocessor and a memory facility in which the equipment's different positions for selecting and fetching samples are stored.

Simply connect a cable and it is ready to go, controlled from the analyser via a serial port (this port is marked "Tecator interface" at the rear of the Sampler). This can be done from any Windows PC. In this way the 5027 Sampler may be integrated with many analytical systems as a reliable and easy-to-use sample presentation machine.

The Sampler enables entire sample batches to be run automatically. The Sampler can handle different sizes of sample cups and is designed to handle long continuous batch runs. As the Sampler is fully automatic, it is extremely easy to use. All the user needs to do manually is to load the sample trays and calibration cups. The Sampler is a durable product and also withstands external environment effects. It operates quietly and moves smoothly.

2.2 Features

- Choice of sample cup sizes: 4 ml, 12 ml, 30 ml or 100 ml
- Turntables for 32, 56, 64 or 120 samples
- Built in safety features
- Rugged and reliable design
- Combines with a wide range of analysis systems
- Very simple to operate
- TTL and RS 232 modes included

2.3 Technical Specification

Power supply:	100-230 V \pm 10%, 50/60 Hz
Power consumption:	20 W
Weight:	14 kg
Dimensions:	450 x 400 x 260 mm
Maximum height with sample probe raised:	370 mm

- Microprocessor controlled unit with all sample positions and operations preprogrammed.
- Features both TTL and RS 232 operations mode.
- Using TTL mode a selector switch determines sample tray type: 32 or 120 positions.
- In the RS 232 mode all sample positions may be addressed randomly and sample probe and turntable movements are controlled from the PC.

3 Unpacking and Checking

Unpack the Sampler with care and check that you have received all parts. If any parts are damaged or missing, please contact your FOSS representative.


Fig. 3:1 1000 5829 Sampler unit with accessory kit (5027-011)

1	Basic unit	5	Power cord, US/Canada type
2	Sample probe	6	Fuse T 0.315 A and Fuse T 0.5 A.
3	Wash bottle	7	User Manual
4	Power cord European type	8	Silicone tube

Depending on choice of sample racks you will receive one of the following sample rack kits:



Fig. 3:2 5000 3697 Accessory kit, 32 cups, 100 ml


1	Sample racks 4 pcs
2	Sample cups, 100 ml 32 pcs

5027044a

Fig. 3:3 1000 9186 Accessory kit, rack cpl with 64 cups, 30 ml

- 1 Sample racks 1 pcs
- 2 Sample cups, 30 ml 64 pcs

5027022a

Fig. 3:4 5000 3713 Accessory kit, rack cpl with 120 cups, 12 ml

- 1 Sample racks 4 pcs
- 2 Sample cups, 12 ml 120 pcs

4 Installation

4.1 General

The 5027 Sampler Unit is designed to be used in a wide range of FOSS Analytical systems. It combines directly with either the Aquatec®, FIAstar®, SCABA, Milko-Scan or WineScan the moment it arrives on your lab bench.

Choose from four different turntables: 32, 56, 64 or a 120 position turntable, and four cup sizes: 4, 12, 30 or 100 ml. Choose the best solution for your lab or use different sample tray configurations for different tasks in the lab. All you have to do is to exchange the sample racks and turn a switch on the control panel. The microprocessor then automatically selects the appropriate function mode from preprogrammed instructions.

4.2 Voltage Feed

Mains Voltage Selector


A mains voltage selector is located at the rear of the sampler and is used for setting the correct mains voltage. You are setting the mains voltage selector with help of a screwdriver or a coin. The selector has the following positions:

Neutral, 100V, 115V, 130V, 200V and 230V, refer to Fig. 4:1.

When the sampler is delivered, the selector is in the neutral position (marked with a line) and there are no fuses fitted. The neutral position means that the primary voltage loop in the sampler is open.

Note:

Never connect the sampler to the mains before setting the mains voltage selector into the desired position and fitting the fuses.

5027026a

Fig. 4:1 Mains voltage selector

Fuses

Two types of fuses are delivered with the sampler. Choose the appropriate fuse in accordance with the following table:

Mains voltage	Fuse
100-130 V	0.5 A
200-230 V	0.315 A

Fit the fuses as shown in Fig. 4:2.

At the rear of the sampler there is the same table as shown on the other page.


Fig. 4:2 How to fit the fuses

4.3 Connecting the Sample Probe

Fit the sample probe by pushing it down into the probe arm. The probe may protrude downwards a maximum of 125 mm.

Fix the probe by tightening the screw on the arm.

For some analysis system different sample probes are used. These are supplied in accessory kits and are described in respective analysis system's User Manual or separate Operator Instructions supplied with the kits.

When the sample probe is fitted, check that the probe arm has slotted into its groove by carefully turning it both to the left and to the right. You should hear a clicking sound to indicate that the arm is in the right position.

For FIAstar systems the Random Access kit needs to be installed, see the enclosed Operators instruction for mounting instructions.

4.4 Fitting the Sample Racks

Fit the sample racks so that the coding on the racks and the centre circle on the sampler match each other.

The sample racks must be positioned correctly to ensure that the numbering of the sample cups will be correct. There will be a space between racks 1 and 4.

4.5 Environmental Conditions

The equipment is designed to be safe under the following conditions:

- Indoor use
- Altitude up to 2000 m.
- Temperature 5 °C to 40 °C.
- Maximum relative humidity 80 % for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C.
- Mains supply voltage fluctuations not exceeding ± 10 % of the rated voltage.
- Transient overvoltage is according to category II, which is normal for this type of equipment.
- Pollution degree 2.

5 Operation

The sampler can be used together with various analysis systems. This section describes how the various systems are used together with the sampler.

Warning

The responsible body shall be made aware that, if the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Modifications, alterations, rebuilding or use of safety parts not authorized by FOSS Analytical AB violates the warranty. FOSS Analytical AB has no responsibility for damages, material or personal, occurring as a result of such actions.

5.1 Connecting Analysis Systems

The connector cable with which you will connect the sampler to the analyzer is usually delivered together with the analysis system.

5.2 Identification of Operating Controls


Fig. 5:1 Operating controls

1	Mains switch	4	Reset button
2	Mains select	5	Rocker switch (functions only in mode 4 and 5)
3	Mode switch		

5.3 Setting the Analysis Configuration


On the rear of the sampler there is a switch marked “Mode”. The “Mode” switch is a hexadecimal switch with which you can select various functions for the sampler, such as different numbers of sample cups or set it up to different analysis systems, see Fig. 5:2.

The “F” switch is a neutral position to be used when the sampler is to be transported or inspected by service personnel. You can change the position of the switch using a small screwdriver. Use the table below to select the appropriate mode of operation:

Mode	Configuration	Trays	Logic
0	Scaba 5600	32	TTL
3	MilkoScan FT 120	56	RS232
4	FIAstar 5010 and Aquatec 5400	32	TTL
5	FIAstar 5010 and Aquatec 5400	120	TTL
9	WineScan and FIAstar 5000	64	RS232
C	FIAstar 5012, Scaba 5610 and 5611	32	RS232
D	FIAstar 5000, FIAstar 5012 and Wine-scan	120	RS232
F	Test and adjustment		TTL

Note:

If you change the position of the switch, the sampler must be reset.

5027028a

Fig. 5:2 The "MODE" switch

5.4 Starting the Sampler

Connect the mains cable to the wall socket and start the sampler by pressing the mains switch on the rear. The green LED at the front lights up when the power is on. The sample probe will automatically move to the wash bottle if it accidentally moved out of this position.

5.5 Operation

Use the table in section 5.3 Setting the Analysis Configuration, to set up the mode switch in the desired position. Remember to restart the sampler for a change to take effect. Usually the PC software of the different analysis systems include a functional or diagnostic test of the sampler, if you wish to check it further.

6 Maintenance

The Sampler does not require any special maintenance apart from the fact that it must always be kept clean.

Note:

Modifications, alterations, rebuilding or use of safety parts not authorized by FOSS Analytical AB violates the warranty. FOSS Analytical AB has no responsibility for damages, material or personal, occurring as a result of such actions.

6.1 Instructions for Cleaning

Clean unit and sample racks using a mild detergent.

7 Accessories and Consumables

Part No	Accessories
1000 5224	Random access kit-5027 Sampler 1 ea Rack for standards 1 ea Communication cable for PC, RS 232, 25-9 pin 10 ea Calibrant cups 1 ea Installation instructions 1 ea Assembly kit
5000 2779	Wash bottle
1000 7361	Communication Cable H9/S25 for PC, RS 232
5000 3573	Sample rack no. 1, 4 ml
5000 3574	Sample rack no. 2, 4 ml
5000 3575	Sample rack no. 3, 4 ml
5000 3576	Sample rack no. 4, 4 ml
5000 3715	Sample rack no. 1, 12 ml
5000 3716	Sample rack no. 2, 12 ml
5000 3717	Sample rack no. 3, 12 ml
5000 3718	Sample rack no. 4, 12 ml
5000 3713	Rack cpl with 120 ml cups, 12 ml
1000 9186	Rack cpl with 64 cups, 30 ml
5000 2996	Sample rack no. 1, 100 ml
5000 2997	Sample rack no. 2, 100 ml
5000 2998	Sample rack no. 3, 100 ml
5000 2999	Sample rack no. 4, 100 ml
5000 3651	Sample probe, 1 channel
1001 0182	Sample probe, 2 channels
1001 0184	Sample probe, 3 channels
1001 3755	Sample probe with filter
1001 3709	Sample probe backflush cpl
1521 0091	Backflush syringe 10 ml
1561 0433	Stopper
5582 0062	Tubing silicone 1/3 mm
1001 3597	Retainer plate
1001 3598	Clamp
1532 0139	Screw M4x30

Part No	Consumables
5522 0034	Sample cups 4 ml 100/pkg
5000 3725	Sample cups 4 ml 1000/pkg
D0399436	Sample cups 30 ml WITH CAPS 750/pkg
5000 3724	Sample cups 12 ml, 2000/pkg
5000 1822	Sample cups 100 ml, 400/pkg
1521 0066	Calibrant bottle with cap
1001 3805	Perforated caps for 12 ml cups, set of 100
1577 0191	Perforated caps for 30 ml cups

8 SAFETY INFORMATION (GB)

8.1 Equipment Ratings

Warning

This device is equipped with a grounding/earthing type power plug for your protection against electrical shock hazard and should only be attached to a properly grounded/earthed receptacle.

Warning

Before replacing the fuses, disconnect incoming mains supply.

Warning

Electrical shock hazard. Covers or panel(s) should be removed by qualified personnel only.

Note:

To maintain the limits for CE approval, only CE approved products may be connected.

8.1.1 Mains Supply

Main unit 100 V/115 V/130 V/200 V/230 V
50-60 Hz

8.1.2 Insulating of External Circuits

The insulation of external, inaccessible circuits is reinforced.

8.1.3 Input and Output Connections

At the rear:

Pulse input: A standard TTL-input.

Tecator interface: A standard RS-232 interface.

8.1.4 Environmental Conditions

The equipment is designed to be safe under the following conditions:

- Indoor use
- Altitude up to 2000 m.
- Temperature 5 °C to 40 °C.
- Maximum relative humidity 80 % for temperatures up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C.
- Mains supply voltage fluctuations not exceeding ± 10 % of the rated voltage.
- Transient overvoltage is according to category II, which is normal for this type of equipment.
- Pollution degree 2.

8.2 Installation

8.2.1 Voltage Feed

Mains Voltage Selector

A mains voltage selector is located at the rear of the Sampler and is used for setting the correct

mains voltage. The mains voltage selector has the following positions: Neutral, 100V, 115V, 130V, 200V and 230V. Refer to the figure.

When the Sampler is delivered, the selector is in the neutral position and there are no fuses fitted.

Note:

Never connect the Sampler to the mains before setting the mains voltage selector into the desired position and fitting the fuses.

Fuses

Choose the appropriate fuse in accordance with the following table:

Mains voltage	Fuse
100-130 V	T 0.5 A
200-230 V	T 0.315 A

Fit the fuses as shown in the figure.


Fig. 8:1 How to fit the fuses

8.2.2 Connecting the Sample Probe

Fit the sample probe by pushing it down into the probe arm. The probe may protrude downwards a maximum of 125 mm.

Fix the probe by tightening the screw on the arm.

For some analysis system different sample probes are used. These are supplied in accessory kits and are described in respective analysis system's User Manual or separate Operator Instructions supplied with the kits.

When the sample probe is fitted, check that the probe arm has slotted into its groove by carefully turning it both to the left and to the right. You should hear a clicking sound to indicate that the arm is in the right position.

For FIAstar systems the Random Access kit needs to be installed, see the enclosed Operators instruction for mounting instructions.

8.2.3 Fitting the sample rack

Fit the sample racks so that the coding on the racks and the centre circle on the sampler match each other.

The sample racks must be positioned correctly to ensure that the numbering of the sample cups will

be correct. There will be a space between racks 1 and 4.

8.3 Operation

The sampler can be used together with various analysis systems. This section describes how the various systems are used together with the sampler.

Warning

The responsible body shall be made aware that, if the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Modifications, alterations, rebuilding or use of safety parts not authorized by FOSS Analytical AB violates the warranty. FOSS Analytical AB has no responsibility for damages, material or personal, occurring as a result of such actions.

8.3.1 Connecting Analysis Systems

The connector cable with which you will connect the sampler to the analyzer is usually delivered together with the analysis system.

8.3.2 Identification of Operating Controls


Fig. 8:1 Operating controls

- | | |
|---|--|
| 1 | Mains switch |
| 2 | Mains select |
| 3 | Mode switch |
| 4 | Reset button |
| 5 | Rocker switch (Functions only in mode 4 and 5) |

8.3.3 Setting the Analysis Configuration


On the rear of the sampler there is a switch marked "Mode". The "Mode" switch is a hexadecimal switch with which you can select various functions for the sampler, such as different numbers of sample cups or set it up to different analysis systems, see Fig. 8:2.

The "F" switch is a neutral position to be used when the sampler is to be transported or inspected by service personnel. You can change the position of the switch using a small screwdriver. Use the table below to select the appropriate mode of operation:

Mode	Configuration	Trays	Logic
0	Scaba 5600	32	TTL
3	MilkoScan FT 120	56	RS232
4	FI Astar 5010 and Aquatec 5400	32	TTL
5	FI Astar 5010 and Aquatec 5400	120	TTL
9	WineScan and FI Astar 5000	64	RS232
C	FI Astar 5012, Scaba 5610 and 5611	32	RS232
D	FI Astar 5000, FI Astar 5012 and Wine-scan	120	RS232
F	Test and adjustment		TTL

Note:

If you change the position of the switch, the sampler must be reset.

5027028a

Fig. 8:2 The "MODE" switch

8.3.4 Starting the Sampler

Connect the mains cable to the wall socket and start the sampler by pressing the mains switch on the rear. The green LED at the front lights up when the power is on. The sample probe will automatically move to the wash bottle if it accidentally moved out of this position.

8.3.5 Operation

Use the table in section 8.3.3 Setting the Analysis Configuration, to set up the mode switch in the desired position. Remember to restart the sampler for a change to take effect. Usually the PC software of the different analysis systems include a functional or diagnostic test of the sampler, if you wish to check it further.

8.4 Instructions for cleaning

Clean unit and sample racks using a mild detergent.

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