IKA® Big Squid

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OPERATING INSTRUCTIONS GB 4

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D

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 61010-1; EN 50081; EN 50082; EN 55011 and EN 60555.

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Staufen, Februar 2004

Reiner Dietsche
Vice President

Wolfgang Buchmann
Corporate Director - Quality Assurance

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Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Contens

	Page
Garantie	2
Safety instructions	4
Correct use	5
Unpacking	5
Commissioning	5
Maintenance and cleaning	5
Technical data	5

Safety instructions

The instrument may not be operated in rooms with explosion hazards. Furthermore, it has to be placed upon a fire-proof and/or not burnable surface.

Attention - Magnetism! Effects of the magnetic field have to be taken into account (e.g. data carriers cardiac pavemakers...).

When using PTFE-coated magnetic bars, the following has to be noted: *Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical systems a temperatures above 300-400°C. Only elementary metals do attack it; halogen hydrocarbons have a reversibly swelling effect.*

Source: Römpps Chemie-Lexikon and „Ullmann“ Vol. 19

Attention: Surface should not be in contact with aggressive medium! Do not place under water!

Correct use

The IKA BigSquid is a magnetic stirrer suitable for mixing thin fluid media up to 1500 ml by using a magnetic bar up to a length of 30mm. The speed can be adjusted between 0 and 1500 1/min.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers:
an IKA BigSquid and operating instructions.

Commissioning

Plug the instrument into a socket with a voltage of 100 ... 240 V \pm 10% and a frequency of 50/60Hz.

Put a magnetic bar into a vessel permeable for magnetic field lines, and fill the vessel with the medium to be processed. Please make sure that the rotary knob for speed adjustment is at the left-hand stop, and then put the vessel on the glass hood. Now you can set the requested speed on the rotary knob.

Maintenance and cleaning

The IKA BigSquid is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning. Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

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