
LEICA MS5 • MZ6 • MZ8 • MZ12

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MIKROSKOP-SYSTEM AB  
Hälsjö 8572  
820 64 Näsviken  
Tel / Fax 0650 - 33211

**Contents**

	<b>Page</b>
<b>Overview</b>	
A conducted tour of the user manual	6
The safety concept must be read by everybody	8
Symbols	10
Controls	13
Short instructions	14
<b>Use</b>	
Setting the interpupillary distance and the viewing height	16
Locating the eyepoint and adjusting the eye contact	17
Changing the magnification	18
Focusing	19
Microscope carrier AX: Axial observation	20
Adjusting the dioptic setting	22
Adjusting for transmitted light	24
Adjusting for incident light	26
Working with the stereomicroscope, tips	28

**Assembly Stands**

Overview: Components of the basic outfit	30
Fitting the column to the incident- and transmitted-light bases	31
Microscope carrier, optics carrier	32
Binocular tubes, optical accessories	33
Transmitted-light stands (bright field)	34
Transmitted-light stands (bright/dark field and EB)	35
Large swinging-arm and table-clamp stands	36
Swinging-arm stand	38

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


## **Illuminators**

6V/10W incident lamp	40
6V/20W incident lamp	42
Coaxial illuminator	44
Fibre-optic light guide	44
Near-vertical illuminator	45
Transformers	46

## **Accessories**

Fitting accessory tubes	48
Phototube: Photography, TV, filming	48
2D- and 3D video systems LEICA 2 and LEICA 3	49
Drawing tube	50
Discussion tube	51
Filter-slide housing	51
Beam splitter: Photography	
Dual-station viewing Cine/TV	52
Attachment for vertical and oblique observation	54
Graticules	54
Stages, polarization	55

## **Notes**

Tips and hints: What to do if ...	56
Care and maintenance of the stereomicroscope, maintenance	58
Calculating total magnification and field diameter	60
Optical data, MS5, MZ6	61
Optical data, MZ8	62
Optical data, MZ12	63

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

## **A conducted tour of the user manual**

### **Do you have a moment?**

We want you to be able to use your stereomicroscope as soon as possible, efficiently and with enjoyment. So you need to get to know the user manual.

### **Safety affects us all**

Before you try to use the instrument, please read this user manual and pay particular attention to the safety notes on pages 8 and 10 and in the individual sections.

### **Are you the one who unpacks and assembles the stereomicroscope?**

On pages 30–55 you will find everything you need to know.

### **Here's how to get to know your stereomicroscope**

Even if you are working with a stereomicroscope for the first time you will soon be able to use it. Study the overview on page 13 to learn what the controls do, and then practice each operation on the basis of the detailed description given in the appropriate section (see contents list on page 4).




The short instructions on page 14 summarize the entire procedure and serve as a checklist.

### **Do you have accessories?**

Stereomicroscopes consist of modules which can be assembled in many different ways. If the configuration of your instrument differs from that of the standard outfit on page 13, refer to the descriptions of the special accessories in the appropriate sections (see contents on page 4). The assembly and operation of the photoautomats are dealt with in a separate user manual.

### **Having problems using the instrument?**

The section on page 56 "Tips and hints: What to do if ..." will help you to find out what you are doing wrong.

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## Safety concept

### General directions

**Before you try to use the instrument, please read the user manual and the safety notes on pages 8 and 10.**

#### • Permitted uses

The LEICA MS5, MZ6, MZ8 and MZ12 stereomicroscopes are optical instruments which use magnification to see the finer details of objects close to the eye. The objects may be either natural, or created by human agency. The basic outfit, which includes a stand and a light source, can be supplemented with various modules for photography, TV, dual-station viewing, and other applications.

#### • Prohibited uses

The use of the instrument in a different manner from that described in this user manual can lead to injury, malfunction and damage.

Do not fit different plugs. Do not dismantle optical systems or mechanical parts unless instructions for doing so are given in the user manual.

#### • Environmental limits to use

The LEICA MS5, MZ6, MZ8 and MZ12 stereomicroscopes are intended mainly for indoor use. If the stereomicroscope is used outdoors, protect it from dirt and damp. Leica lamps, and Leica stands which include lamps, may not be used outdoors.

#### • Servicing

Repairs may only be carried out by Leica-trained service technicians. Only original Leica spare parts may be used.

#### • Responsibilities of person in charge of instrument

Ensure that personnel who use the instrument have read and understood this user manual and in particular the safety instructions on pages 8 and 10.



Ensure that the LEICA MS5, MZ6, MZ8 and MZ12 stereomicroscopes are operated, maintained and repaired only by authorized and trained personnel.

#### • Legal requirements

Adhere to general and local regulations relating to accident prevention and environmental protection.

#### • Conformity with European Community directive

The LEICA MS5, MZ6, MZ8 and MZ12 stereomicroscopes and their accessories are constructed in accordance with the latest technologies and are provided with a statement of conformity with EC requirements.

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## Safety concept

You will find these symbols in the user manual:

### Safety direction

This symbol denotes information which must be read and obeyed.

Failure to adhere to these directions can cause injury.

### Warning: Electrical current

This symbol denotes information which must be read and obeyed.

Failure to adhere to these directions can cause injury.

### Important warning

This symbol denotes information which requires particular attention.

Failure to adhere to these directions can cause malfunctioning in, or damage to, the instrument.

### Important information


This symbol denotes additional information or explanations which promote understanding.

### Action

⇒ This symbol within the text indicates that certain operations must be carried out.

### Explanatory notes

- This symbol within the text stands for additional information and explanations.

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**In this overview**

we will show you the controls and the functions of the standard outfit.

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
## **Stereomicroscope**

### **1 Magnification changer**

Drive knobs with magnification scale.  
– for MS5: 5 steps  
– for MZ6: 6:1 zoom  
– for MZ8: 8:1 zoom  
– for MZ12: 12.5:1 zoom

### **2 Engaging ratchet positions**

– for MZ6: switch engages 7 magnification steps  
– for MZ8: milled ring engages 8 magnification steps  
– for MZ12: switch engages 10 magnification steps

### **3a Focusing drive**

#### **3b Coarse and fine focusing**

Larger knob: coarse focusing  
Smaller knob: fine focusing.

#### **3c Adjustable ring**

Changes the ease of movement of the coarse/fine drive.

#### **4 Securing screw**

Holds the optics carrier in the microscope carrier.

#### **5 Securing screw**

Holds the binocular tube (or accessory) to the optics carrier.

#### **6 Interchangeable objective**

Various objectives (achromatic, plan-achromatic, planapochromatic), depending on the outfit.

## **Binocular tube**

### **7 Adjustable eyepiece tubes:**

Interpupillary distance adjustable within the range 52mm – 76mm.

### **8 Wide-angle eyepieces for spectacle wearers**

Dioptric setting adjustable within the range from +5 to –5, and eyecups adjustable.

### **9 Securing screws**

These hold the eyepieces in the binocular tube.

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### Short instructions Checklist

#### In this overview

we have summarized the working sequence. You will find detailed user instructions on the pages indicated.

#### Inspect the cleanliness

- ⇒ Remove dust and dirt from the instrument (page 59).
- ⇒ Cleaning: Optical accessories, glass stage plates, and any graticules used (page 59).

#### Only for binocular tube with variable angle of observation

Adjust viewing height:

- ⇒ Swing the eyepiece tubes upwards or downwards (page 16).

#### Only for photomicrographic outfits

- ⇒ Engage observation beam path: Lever at position "VIS" (page 48).
- ⇒ For trinocular tube: Pull diaphragm lever out to the left.
- ⇒ Open double-iris diaphragm (position "OPEN").

#### Adjust the interpupillary distance

- ⇒ Using both hands, displace the eyepiece tubes (page 16).

#### Select the degree of eye contact

- ⇒ Secure the eyecup in the position required (page 17).

#### Adjust the dioptic setting

- ⇒ Adjust the dioptic setting precisely for each eyepiece (page 22).

#### Set the magnification required


- ⇒ Turn the rotary knob to the desired magnification (page 18).

#### Focusing

- ⇒ Look into the eyepieces and focus on the object (page 19).
- Only for coarse/fine drive:  
Focus approximately with the large knob, and then accurately with the small knob.
- ⇒ If the focusing movement is too loose or too tight, adjust the ease of movement (page 19).

#### Microscope carrier AX

- ⇒ For measuring, photography and polarization, engage the vertical beam path (page 20).

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## Use

### Viewing height, interpupillary distance, eye contact

#### In this section

you will learn how to set up the microscope so that you can work in a relaxed and comfortable manner, largely free of fatigue.

### Only for the binocular tube with variable angle of observation, 10°–50°

- The viewing height is correct when your head and back are comfortably positioned.
- ⇒ Holding the eyepiece tubes at their bases, swing them up or down (see arrow).

### Adjusting the interpupillary distance

- The interpupillary distance is correct when you see a single circular field when viewing with both eyes.
- ⇒ Look into the eyepieces and use both hands to move the eyepiece tubes closer together or further apart.

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### Locating the eyepoint

- For 10x/21B wide-field eyepieces for spectacle wearers, the distance from eye to eyepiece is 22mm.
  - When you can see the entire field of view without restriction, you have found the correct position (eyepoint).
- ⇒ To find the eyepoint, slowly bring your eyes closer towards the eyepieces.

### Selecting the degree of eye contact

#### If you do not wear spectacles, and want closer contact to the eyepieces

- ⇒ Holding the diopter ring, turn the eyecup anticlockwise until it is released.
- ⇒ Pull the eyecup upwards.
- ⇒ Holding the diopter ring, turn the eyecup clockwise.

### Spectacle wearers

need the whole of the available distance. Move the eyecups to the lowest position.

### Adjusting the dioptic setting and the parfocality

You must adjust the dioptic setting for both eyes exactly as described on page 22. Only then will the image be parfocal, i.e. it will remain sharp whenever the magnification is changed and there will be no need to refocus. Make the most of this opportunity, which is not offered by every stereomicroscope.

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- Las dioptrías se ajustan de +5 a -5.

## Use

### Changing the magnification, setting the working distance, focusing

#### In this section

we will show you how to change the magnification, to engage ratchets at certain magnification steps, and to focus.

#### Magnification changer

MS5: 5 steps  
MZ6: 6:1 zoom  
MZ8: 8:1 zoom  
MZ12: 12.5:1 zoom

- Zoom means continuous change of magnification. Ratchets can be engaged for certain magnifications (see below).

#### Changing the magnification

- ⇒ Look into the eyepieces.
- ⇒ Focus on the object (page 19).
- ⇒ Turn the magnification changer to the position required.
- The magnification changer can be operated from either the right or the left.
- The scale on each drive knob is inscribed with the zoom magnification factors. These apply when a 1.0x objective is used.

#### Ratchet steps

MZ6: 0.8, 1.0, 1.25, 1.6, 2.0, 2.5, 3.2  
MZ8: 0.8, 1, 1.25, 1.6, 2, 2.5, 3.2, 4  
MZ12: 1, 1.25, 1.6, 2, 2.5, 3.2, 4, 5, 6.3, 8

- ⇒ For MZ6 and MZ12, press switch
- ⇒ For MZ8, turn knurled ring anticlockwise.

#### Set the ratchets

- Engage the ratchets whenever magnifications need to be exactly reproducible.

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## Scale of magnifications

Only for MZ8:  
Other adhesive scales can be applied to the rotary knobs for other eyepiece/objective combinations (see page 32).

## Setting the working distance, focusing = obtaining a sharp picture

- ⇒ Place the object beneath the objective.
- ⇒ Engage the lowest magnification (page 18).
- ⇒ Look into the eyepieces.
- ⇒ Using the rotary knob, bring the object into focus.
- To focus, use the focusing drive to raise or lower the stereomicroscope until the feature of interest is in focus. The working distances of the various objectives are listed on page 59.
- The focusing drive can be operated from either the right or the left.

## Only for focusing drive coarse/drive

- ⇒ Focus approximately with the larger knob.
- ⇒ Focus accurately with the smaller knob.

## Adjusting the free running of the focusing drive

Does the focusing drive run too slack or too tight, or does the equipment drift downwards?  
You can match the ease of movement to the weight of the outfit and to your personal requirements.

## Focusing drive

- ⇒ Grip the drive knobs firmly with both hands and turn them in opposite directions until you obtain the required degree of resistance during focusing.

## Focusing drive, coarse/fine

Turn the black ring on the right-hand drive knob until you obtain the required degree of resistance during focusing.

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### **Microscope carrier for stereoscopic observation**

This microscope carrier permits only stereoscopic observation, using both beam paths.

- The optics carrier can be turned sideways in the microscope carrier if the user needs to work from the side:

- ⇒ Release the clamping screw
- ⇒ Turn the optics carrier sideways as required
- ⇒ Retighten the clamping screw.

The microscope carrier can be attached to the focusing drive in either of two positions:

- For short working distances and for flat objects, use Position 1. Use this position also in conjunction with the attachment for vertical and oblique observation (see detailed instructions).
- Use Position 2 for long working distances and for bulky objects.

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*Optics carrier turned  
sideways in microscope  
carrier.*

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## Microscope carrier AX

### Only for MS5 and MZ6

The MS5 and MZ6 stereomicroscopes can be equipped with the microscope carrier AX for stereoscopic and axial observation instead of the standard microscope carrier.

- The optics carrier **cannot** be turned sideways when in the microscope carrier AX (see page 20).
- The microscope carrier AX can only be fitted **in one position** to the focusing drive (see page 20).

## Stereoscopic observation

A three-dimensional image is essential for studying and working on spatial objects.

- ⇒ Shift the stereomicroscope to the central position.

## Axial observation

The parallax-free imaging produces more precise results in photography, measuring and polarization.

- ⇒ Shift the stereomicroscope to the left or the right position. The object is now located directly beneath either the left or the right beam path.
- When photographing or measuring, select the beam path above which the phototube or measuring graticule is located.

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#### In this section

you will learn

- how to adjust the dioptic settings for each eye, and
- how to set up the stereomicroscope parfocally (refer also to page 17).

Each user needs to carry out these operations once.

#### A Preparations

- ⇒ Set up the illumination (pages 24–27).
- ⇒ Set the interpupillary distance and, if necessary, the viewing height (page 16).
- ⇒ Ensure that the diaphragms, including double-iris diaphragm are open (VIS and OPEN on phototube).
- ⇒ Adjust the eyecups (page 17).
- ⇒ Set both eyepieces to zero.
- ⇒ Set the microscope carrier AX for stereoscopic observation.

#### B Focus on the test object

##### 1. Flat test object

- ⇒ Place such an object beneath the objective.

##### 2. Lowest magnification

- ⇒ Engage the lowest magnification position (page 18).
- ⇒ Use the coarse focusing drive to set the working distance (page 19).

##### 3. Look into the eyepieces

- ⇒ Using the focusing drive, bring the test object into focus.

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##### 5. Image sharpness

- ⇒ Using the focusing drive, bring the image into sharp focus.

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### C Adjusting the dioptric settings

#### 6. Lowest magnification

- ⇒ Engage the lowest magnification position.

#### 7. Do not look into the eyepieces

- ⇒ Rotate the eyepieces anticlockwise in the "+" direction as far as the stop (+5 diopters).

#### 8. Look into the eyepieces

- ⇒ Slowly rotate each eyepiece individually in the "-" direction until each eye sees the object sharply imaged.
- ⇒ Record the dioptric setting for each eye.

#### 9. Eyepieces

- ⇒ Reset the eyepieces to "0".
- ⇒ Repeat steps 4 to 8 twice more and make a note of the dioptric settings.
- ⇒ Finally, set the median values on the eyepieces.

### D Inspecting the parfocality

#### 10. Highest magnification

- ⇒ Engage the highest magnification position.

#### 11. Refocus

- ⇒ if necessary.

#### 12. Magnification changer

- ⇒ Slowly move the magnification changer from the lowest magnification to the highest.
- The image sharpness must remain constant. Otherwise, repeat the procedure.

### E Graticule (reticle) in eyepiece?

The procedure for adjusting the dioptric settings and the parfocality is slightly different if a graticule is in position. Refer to the user manual for the MPS modular photographic systems or for measurement accessories.

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## Use

### **Illumination: Transmitted light**

Leica transmitted-light stands are not to be used outdoors.

### **Transmitted-light stand for bright and dark field**

Switch on the illuminator at its source (see user manual for illuminator).

### **Transmitted light, bright field**

- ⇒ To obtain bright-field illumination, push the knob fully in.
- Bright field is suitable for examining translucent objects having contrasty structures. The object is directly illuminated from below and is seen in its natural colours against a bright background.

### **Transmitted light, dark field**

- ⇒ To obtain dark-field illumination, pull the knob fully out.
- In the dark-field technique, the light from the source meets the object very obliquely from below. Only rays which are reflected or refracted at interfaces reach the objective, so edges and structural elements appear bright on a dark background, even if the contrast of the original object was extremely poor.

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## Transmitted-light stand bright field, 20W

Follow the directions for using the built-in transformer (page 34).

- ⇒ Using the thumb wheel switch on the rear of the stand, switch on the transformer and adjust the brightness.

## Transmitted- or incident light

- You can also connect an incident lamp to the transmitted-light stand (page 33).
  - Either transmitted light or incident light can be engaged.
- ⇒ For transmitted light: Switch in position I.  
⇒ For incident light: Switch in position II.  
⇒ Use thumb wheel to adjust brightness.

## Transmitted light, bright field or inclined transmitted light

- ⇒ Switch on transmitted light (see above).

## Transmitted light, bright field

- ⇒ For transmitted light, bright field: Pull lever forwards.
- Bright field is suitable for examining translucent objects having contrasty structures. The object is directly illuminated from below and is seen in its natural colours against a bright background.

## Inclined transmitted light

- ⇒ Move the lever slowly towards the column of the stand until the desired effect is obtained.
- Inclined transmitted light which traverses the object obliquely will produce effects advantageous for observing small translucent objects such as foraminifera and fish eggs.

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#### **Incident lamp, 10W**

Follow the directions for using the transformer (page 46).

⇒ Switch on the transformer and use the rotary knob to adjust the brightness (page 46).

⇒ Hold the lamp housing and release the clamping screw.  
⇒ Position the light spot on the object.  
⇒ Tighten the clamping screw.

⇒ Swing in the heat-absorbing filter.

⇒ Enlarge or reduce the light spot by turning the socket of the bulb.

- The 6V/10W lamp is a source of oblique incident light for illuminating spatial objects.
- To reduce or eliminate shadows, use two lamps.
- The colour temperature at nominal voltage (6V) is 2700K.

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## Incident lamp, 20W

Follow the directions for using the transformer (page 46).

- ⇒ Switch on the transformer and use the rotary knob to adjust the brightness.

## Illuminating the object

- ⇒ Push the condenser lens forwards or backwards until the slider engages.
- ⇒ Hold the lamp housing and release the clamping screw.
- ⇒ Position the light spot on the object.
- ⇒ Tighten the clamping screw.
- ⇒ Match the size of the lamp spot to requirements by pushing the condenser lens forwards or backwards.

- ⇒ Swing in a filter if required (page 42).
- The 6V/20W lamp is a source of oblique incident light for illuminating spatial objects.
- To reduce or eliminate shadows, use two lamps.
- The colour temperature of 3200K at nominal voltage (6V) is ideal for photography. Set the regulating transformer to nominal voltage (page 46).

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## Use

### Working with the stereomicroscope

#### In this section


we will show you how to work smoothly with your stereomicroscope and how to study objects routinely.

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
### Preconditions for working comfortably

- Pay great attention when setting up your stereomicroscope. You will only benefit fully from the outstanding optical and environmental advantages of the instrument if you set it up precisely in accordance with the directions already given on pages 16 to 23.
- Arrange your workplace as well as possible, and pay attention to the heights of chair and table.
- When sitting, use the entire seat area and also the back support of the chair.
- Support your forearms.
- During work breaks, stand up and move around for exercise.

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### Tips for your work

The stereomicroscopes are parfocally matched. This means that, once you have focused on a feature in the manner described, taking account of dioptic settings (page 22), you can zoom the magnification up and down without refocusing. You only need to refocus if you want to look at higher or lower planes in the object.

- When focusing on a fresh object, always start with the lowest magnification, because the larger field of view makes the object easier to find. Afterwards, select the magnification you want.

- ⇒ Place the object beneath the objective.
- ⇒ Illuminate the object.
- ⇒ Select the lowest zoom position.
- ⇒ Using the focusing drive, set the working distance approximately (page 19).
- ⇒ Look into the eyepieces and use the focusing drive to bring the object precisely into focus.

- ⇒ Select the highest zoom position.
- ⇒ Use the focusing drive to refocus.

**Now the object should remain in focus at all zoom magnifications. Otherwise, check that your dioptic settings have been adjusted exactly as described on page 22.**

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