## Assembly

## Montage

## Aufbau

## Montaje

### The components, stands

- 1a Incident-light base with stage plate
- 1b Transmitted-light base with glass stage plate
- 2a Side-faced column with drive housing coarse/fine
- 2b Side-faced column with focusing drive
- 3 Microscope carrier
- 4 Optics carrier
- 5 Interchangeable objective
- 6 Binocular tube
- 7 Wide-field eyepieces for spectacle wearers, with eyecups

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
### In the following sections

we will explain how to assemble a complete stereomicroscope outfit, including the incident- and transmitted-light stands. Detailed directions about the transmitted-light stands, the swinging-arm stand and the fitting of the remaining accessories are all described in the appropriate sections (see contents on page 4).

#### Side-faced column → incident-light base

- ⇒ Remove the stage plate.
- ⇒ Introduce 3 Allen screws into the underside of the baseplate and secure them to the bottom end of the side-faced column.
- ⇒ Insert the stage plate.

To fit the microscope carrier and the stereomicroscope, see page 32.

#### Side-faced column → transmitted-light base

**Transmitted-light stand, bright field, 20W: Pull the power cable out of the socket in the stand.**

- **Transmitted-light stand for bright and dark field:** Pull the fibre-optic light guide out of the connector in the stand.
  - ⇒ Remove the glass stage plate.
  - ⇒ Carefully invert the transmitted-light base and remove 7 Allen screws from the base.
  - ⇒ Lift off the lower part of the base.
  - ⇒ Reverse the upper part of the base.
  - ⇒ Introduce 3 Allen screws into the cover plate from beneath and secure them to the bottom end of the side-faced column.
  - ⇒ Replace and secure the base.
  - ⇒ Insert the glass stage plate.
  - ⇒ Fit the microscope carrier and stereomicroscope as described on page 32.

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## Assembly

### Microscope carrier, optics carrier

#### In this section

we will explain how to fit the microscope carrier and the optics carrier to the incident- and transmitted-light stands. The combination of the optics carrier with the inclinable focusing drive is described on page 39.

### Microscope carrier → focusing drive

The microscope carrier can be fitted in either of two positions to the focusing drive:

- ⇒ Use an Allen screw to secure the microscope carrier in the position required.

### Microscope carrier AX

- The microscope carrier AX can be fitted to the focusing drive only in position 1 (see above).
- The optics carrier cannot be turned sideways in the microscope carrier AX.

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- ⇒ Gently lower the optics carrier into the microscope carrier.
- ⇒ Use the screw to secure the optics carrier in the position required.

### Apply the adhesive scale of magnifications

- Only for MZ8:
- ⇒ Remove the screws.
  - ⇒ Remove the rotary knobs.
  - ⇒ Apply the adhesive scales.
  - ⇒ Replace the rotary knobs.

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## Objective

⇒ Screw on the selected objective anticlockwise.

### Only for MZ8:

- The 1.0x plano objective has a standard thread accepting commercially-available photographic accessories.
- To use the 1.0x plano objective with the coaxial illuminator, with the near-vertical illuminator, or with the attachment for vertical and oblique observation, fit one or two spacer rings to the optics carrier.
- If you do not want to use such accessories as the analyzer (page 55) or the fibre-optic ring lamp, then for some objectives you can unscrew and remove the spacer ring to improve access to the object:
  - for the MS5, MZ6, MZ8: 1.0x Planapo and 1.6x Plan
  - for the MZ12: 1.6x Planapo

## Binocular tube

- ⇒ Release the clamping screw.
- If you have an accessory such as a coaxial illuminator (page 44), a phototube (page 48) or a filter-slide holder (page 51), fit it now to the optics carrier. Tighten the clamping screw.
- ⇒ Tilt the binocular tube and insert the edge of its dovetail ring beneath the two lugs on the optics carrier.
- ⇒ Move the binocular-tube from side to side until the locating screw on the optics carrier engages the guide.
- ⇒ Tighten the clamping screw.

## Eyepieces

- ⇒ Insert a graticule if appropriate (page 54).
- ⇒ Release the clamping screws on the eyepieces.
- ⇒ Insert the eyepieces **fully** into the eyepiece tubes. Check that they have seated correctly.
- ⇒ Tighten the clamping screws.
- You can extend the overall magnification range by using 16x, 25x or 40x wide-field eyepieces for spectacle wearers.

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## Assembly

### Transmitted-light stands

#### Transmitted-light stand, bright field, 20W

- ⇒ Attach the side-faced column to the transmitted-light base (page 31).
- ⇒ Fit the microscope carrier and the stereomicroscope (page 32).

#### Fitting the bulb

- Do not touch the 6V/20W halogen bulb directly with your fingers, because fingermarks can cause it to crack when heated. Remove fingermarks immediately with a cloth dipped in pure alcohol.
- ⇒ Extract the slider together with the bulb socket.
- ⇒ Hold the new bulb with a cloth and insert it fully into the mount.
- ⇒ Push the slider carefully but completely back along the rail to the stop.

#### Filters

- The transmitted-light stand contains a KG1 filter.
- You can insert other 50mm diameter filters if you need them:
- ⇒ Pull out the slider.
- ⇒ Insert the 50mm diameter filter.
- ⇒ Push the slider carefully but completely back along the rail to the stop.

#### Starting up

- ⇒ Set the available voltage on the voltage selector:
  - Position 115 for 100V–120V.
  - Position 230 for 200V–240V.
- ⇒ Connect the power cable to the socket of the base and to the power source.
- To switch on the transformer and operate the transmitted-light stand, see page 25.

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### Changing fuses

- ⇒ Disconnect the power cable from the stand.
- ⇒ Pull out the fuse holder.
- ⇒ Fuse at rear: Press out downwards.
- ⇒ Spare fuse at front: Pull out sideways and insert into holder at rear.
- ⇒ Push the fuse holder back into position again.

### Incident lamp

- ⇒ Connect an incident lamp if required.
- To switch on and off, see page 25.

### Transmitted-light stand for bright and dark field

- ⇒ Fit the side-faced column to the transmitted-light base as described on page 31.
- ⇒ Fit the microscope carrier and the stereomicroscope as described on page 32.
- ⇒ Connect the Volpi fibre-optic light guide to the stand and to the light source.
- To switch on and operate, see page 24.

### Starting up

Before using the light source, read the user manual for it.

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- ⇒ Introduce 3 Allen screws into the adapter from beneath and secure it to the bottom end of the side-faced column.
- ⇒ Use 2 Allen screws to secure the adapter and side-faced column together to the transmitted-light base EB.

**Detailed information** about starting-up and use is given in the separate user manual for the transmitted-light base EB.

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## Setting up

### Large swinging-arm and table-clamp stands

- ⇒ **Table-clamp stand:** Use the clamp to secure the column to a table top 20mm – 50mm thick.
- ⇒ **Swinging-arm stand:** Introduce 4 Allen screws into the base from beneath and secure them to the end of the column.

### Safety ring → column

The safety ring secures the instrument on the column and must always be positioned beneath the horizontal arm.

- ⇒ Position the safety ring with the screw head uppermost and towards you, and slide it down the column.
- ⇒ Tighten the clamping screw.

### Horizontal arm → column

- ⇒ Lower the horizontal arm on to the safety ring, with its end pointing towards you.
- ⇒ Tighten the clamping screw.

### Limitation

The screw head must engage the notch. It limits the lateral movement. It may not be removed, and neither may the safety ring be milled down. If the horizontal arm is swung out beyond the limit, the outfit may fall over.

### Drive housing → horizontal arm

Here follow directions for fitting the drive with side-faced carrier rod. To fit the drive housing with peg, refer to page 38.

- ⇒ Release the other clamping lever.
- ⇒ Tilt the carrier rod downwards.
- ⇒ Retighten the clamping lever.

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- ⇒ Slide the drive housing with side-faced carrier rod on to the horizontal arm.
- ⇒ Tighten the clamping lever.

The clamping levers can be positioned as required:

- ⇒ Pull out the clamping lever along its axis and turn it.

#### Microscope carrier → drive housing

- ⇒ Secure as described on page 32.

#### Optics carrier → microscope carrier

- ⇒ Secure as described on page 32.

#### Possibilities for adjustment

##### Tilting

- ⇒ Release the clamping lever.
- ⇒ Tilt the instrument.
- ⇒ Retighten the clamping lever.

##### Focusing

- ⇒ Set the distance approximately with the larger knob.
- ⇒ Fine-focus with the smaller (concentric) knob.
- ⇒ Using the ring on the right-hand drive knob, adjust the ease of movement of the drive.

#### Changing the working distance

- ⇒ Release the wing screw.
- ⇒ Raise or lower the horizontal arm.
- ⇒ Retighten the wing screw.

Position the safety ring (see page 36)

#### Swinging out sideways

- ⇒ Release the wing screw.
- ⇒ Displace the horizontal arm.
- ⇒ Retighten the wing screw.

When swinging sideways, respect the end-stops (see page 36)

#### Changing the overhang

- ⇒ Release the wing screw
- ⇒ Pull out or push in the horizontal arm
- ⇒ Tighten the wing screw.

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Only for MS5, MZ6, MZ8

#### Assembly

- ⇒ Introduce 4 Allen screws into the base from beneath and secure them to the bottom end of the column.

#### Safety ring → column

The safety ring secures the instrument and must always be fixed beneath the cross-member.

- ⇒ Slide the safety ring down the column.
- ⇒ Tighten its clamping screw.

#### Cross-member → column

- ⇒ Lower the cross-member on to the safety ring.
  - ⇒ Tighten the clamping lever.
- The clamping levers can be positioned as required:
- ⇒ Pull out the clamping lever along its axis and turn it.

#### Horizontal arm → cross member

- ⇒ Slide the horizontal arm into the cross-member with the connecting piece towards you.
- ⇒ Tighten the clamping lever.

#### Drive housing → horizontal arm

- ⇒ Insert the peg of the inclinable focusing drive into the connecting piece.
- ⇒ Tighten the wing screw.

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- The drive housing / microscope carrier can be fitted to the horizontal arm in either of two positions.
- The round port in the microscope carrier is designed to accept a light source.

### Optics carrier → microscope carrier

- ⇒ Fit the optics carrier in the microscope carrier.
- ⇒ Turn the optics carrier sideways in the microscope carrier as required.
- ⇒ Tighten the clamping screw.

### Possibilities for adjustment

#### Tilting

- ⇒ Release the Allen screw.
- ⇒ Tilt the instrument.
- ⇒ Retighten the Allen screw.

#### Focusing

The focusing knobs are located on the inclinable drive housing.

Is the focus movement too light or too heavy? Does the outfit tend to slide downwards? You can adjust the movement in accordance with the weight of the outfit and with your own personal requirements:

- ⇒ Hold the drive knobs with both hands and turn them in opposite directions until the ease of movement is satisfactory.

#### Altering the working distance

- ⇒ Release the clamping lever.
- ⇒ Raise or lower the cross-member as required together with the horizontal arm.
- ⇒ Retighten the clamping lever.

#### Reposition the safety ring (see above)

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### 6V/10W incident lamp

#### Fitting

- ⇒ Release the clamping screw.
- ⇒ Unfold the lampholder arm.
- ⇒ Tighten the clamping screw.

#### Lampholder → clamp for side-faced columns

- ⇒ Use the screw to secure the clamp to the side-faced column, either above or below the drive housing.
- ⇒ Use the connector to attach the lampholder arm to the clamp, on either the right or the left.

#### Lamp housing → lampholder arm

- ⇒ Attach the lamp to the connector on the lampholder arm.

- ⇒ Fit the filterholder, complete with KG1 heat-absorbing filter, to one of the pegs.
- A total of 4 filterholders with 50mm diameter filters can be attached.

### Operating the lamp

#### Follow the directions for using the transformer (page 46).

- ⇒ Connect the cable between the lamp housing and the transformer (page 46).
- ⇒ Switch on the transformer (page 47).
- To operate the lamp, see page 26.

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## Changing the bulb

**Disconnect the lamp cable from the transformer.**

- You can separate the lamp housing and the bulb socket. Line up the arrow head on the lamp mount with its shaft on the lamp housing.
  - Do not touch the 6V/10W halogen bulb directly with your fingers, because fingerprints can cause it to crack when heated. Remove fingerprints immediately with a cloth dipped in pure alcohol.
- ⇒ Turn the bulb socket to the separation position.
- ⇒ Pull the bulb socket forcibly out of the housing.
- ⇒ Hold the new bulb in a cloth and push it into the socket.
- ⇒ Replace the socket and bulb together in the housing. The socket clicks into position.

## Other possibilities for fitting

### Lampholder → cast foot or transmitted-light stand EB

- ⇒ Fit the 25mm diam. adapter to the cast foot or to one of the short posts on the base of the transmitted-light stand EB.
- ⇒ Using the connector, attach the lampholder to the adapter.

### Lampholder → objective

- ⇒ Remove the objective.
- ⇒ Slide the adapter over the objective mount.
- ⇒ Using the connector, attach the lampholder to the adapter.

### Lampholder → inclinable drive housing

- ⇒ Unscrew a connector from the lampholder.
- ⇒ Release the clamping screw.
- ⇒ Attach the thread of the lampholder directly to the drive housing.
- ⇒ Retighten the clamping screw.

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### 6V/20W incident lamp

#### Fitting

- ⇒ Release the clamping screw.
- ⇒ Unfold the lampholder arm.

- ⇒ Secure the lampholder to the lamp.

- ⇒ Use the screw to secure the clamp to the side-faced column, either above or below the drive housing.
- ⇒ Use the thread connector to attach the lampholder arm to the clamp, on either the right or the left.
- The thread connector also enables the lampholder to be attached to the cast foot or to one of the short posts on the EB base.

#### Adjusting the lamp

- ⇒ Position the lamp
- ⇒ Tighten the clamping screw.
- ⇒ Push the filterholder, complete with diffusing filter, on to one of the pegs, and swing it in.
- A total of 4 filterholders with 32mm diameter filters can be attached. A heat-absorbing filter is built in.

#### Operating the lamp

Follow the directions for using the transformer (page 46).

- ⇒ Connect the cable between the lamp housing and the transformer (page 46).
- ⇒ Switch on the transformer (page 47).
- To operate the lamp, see page 27.

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- Do not touch the 6V/20W halogen bulb directly with your fingers, because fingermarks can cause it to crack when heated. Remove fingermarks immediately with a cloth dipped in pure alcohol.
- ⇒ Press the spring clips together and pull off the rear part of the housing.
- ⇒ Hold the new bulb in a cloth and push it into position.
- ⇒ Press the spring clip together and close the housing again.
- ⇒ Centre the bulb (see below).

### Centring the bulb

The 6V/20W lamp is a powerful source of illumination. The more carefully you set it up, the more uniform will be the lighting. For this reason, after fitting a new bulb, recentre it.

- ⇒ Switch on the lamp.
- ⇒ Swing out the diffusing filter.
- ⇒ Pull the condenser lens fully out.
- ⇒ Hold the lamp and release the clamping screw.
- ⇒ Point the lamp at a sheet of paper so that two filament images are seen.
- ⇒ Tighten the clamping screw.
- ⇒ Using the slotted screw, bring the two filament images into focus.
- ⇒ Displace the filament images so that they are opposite to one another, but in contact.
- ⇒ Swing in the diffusing filter once more.

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**Coaxial illuminator**

- This illuminator is for observing and photographing flat, highly-reflecting objects such as wafers, integrated circuits, and polished metal sections.
  - Please refer to the detailed directions in the user manual M2-232-004.
  - The coaxial incident light housing fits like an accessory tube between the optics carrier and the binocular tube, or between the optics carrier and the phototube (page 48).
  - Only for MZ8: If you use a spacer ring with the 1.0x plano objective you will be able to illuminate the largest field of view completely.
  - A quarter-wave plate is required when the microscope carrier AX is used in the vertical position.
- ⇒ Using the clamping screw, secure the quarter-wave plate to the objective mount.
  - ⇒ Using the milled ring, turn the quarter-wave plate until the required effect is obtained.

**Fibre-optic light guide**





- ⇒ Unscrew the clamping screw from the clamp.
- ⇒ Hold the clamp in the port of the microscope carrier.
- ⇒ Push the clamping screw through the holes in the microscope carrier and in the broader part of the clamp, and secure it finger-tight to the narrower part of the clamp.
- ⇒ Place the light guide in the clamp.
- ⇒ Switch on the light source.
- ⇒ Position the light spot on the object.
- ⇒ Tighten the clamping screw.
- The clamp for fibre-optic light guides cannot be attached to the microscope carrier AX.

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## Near-vertical illuminator

Only for MS5, MZ6, MZ8:

- The near-vertical (approx. 5°) light beam illuminates depressions and cavities.
- Shadowing by tools is eliminated.
- This illuminator can only be used with the achromatic objectives having focal lengths 100mm, 150mm, 175mm and 200mm.

- ⇒ Screw off the objective.
- ⇒ Release the clamping screw and remove the adapter ring from the incident-light housing for near-vertical light.

For MS5, MZ6:

- ⇒ Screw the adapter ring on to the microscope carrier.

For MZ8:

- ⇒ First screw two spacer rings on to the microscope carrier, then screw the adapter ring to the lower of the two spacer rings.
- ⇒ Secure the incident-light housing to the adapter ring.
- ⇒ Screw the objective into position (see above).
- ⇒ Release the left or right clamping screw and change the cover plate accordingly.
- ⇒ Attach the fibre-optic light guide on either the left or the right.

To switch on the illuminator, refer to the user manual provided with it.

### No light?

- ⇒ Turn the rotary knob.
- A movable prism directs the light from the left or right outlet on to the object.

- Clamp with fibre-optic light guide
- Near-vertical illuminator
- Bride avec guide-lumière à fibres optique
- Eclairage vertical
- Klemme mit Glasfaser-Lichtleiter
- Vertikalbeleuchtung
- Pinza con conductor de luz de fibra óptica
- Iluminación vertical

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
- 1 Voltage selector
- 2 Connections for two lamps (maximum combined loading 40W)
- 3 Connection for power cable
- 4 Fuse holder F1.0A
- 5 Power switch  
○ off  
| on
- 6 ON display: Instrument is switched on
- 7 Regulating knob
- 8 6V display: The maximum voltage has been reached.

### Operation

⇒ Set the available voltage on the voltage selector:

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

- ⇒ Connect the power cable to the socket and to the power supply.
- ⇒ Connect the lamp. The sockets accept two lamps with a maximum combined loading of 40W.
- ⇒ Turn on the power switch. The diode shines.
- ⇒ Use the rotary knob to adjust the brightness. Above 6V the diode shines.

### Switching on and off

You can prolong the life of the bulb if you:

- ⇒ first switch off the lamp by turning the rotary knob into the black range
- ⇒ and only then switch the transformer on or off.

### Changing fuses

Fuses 1.0A, 5x20 (stock no. 167 651)

- ⇒ Lamp off.
- ⇒ Power switch off.
- ⇒ Turn the fuse holder anticlockwise and extract it.
- ⇒ Fit a new fuse and replace the fuse holder.

### Overload

- ⇒ Lamp off.
- ⇒ Power switch off.
- ⇒ Wait a few seconds.
- ⇒ Power switch on.
- ⇒ Lamp on.

### Step transformer

- You can attach a 6V/10W lamp to the step transformer.
- ⇒ Set the voltage selector on the rear side to 115V or 230V (page 46).
- ⇒ Attach the power cable to the power supply.
- ⇒ Connect the lamp.
- ⇒ Adjust the brightness with the three-step switch (4V/5V/6V).

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#### Phototubes, TV, video, photography

##### Fitting

- All accessory tubes which fit between the optics carrier and the binocular tube are fitted in the same manner.
  - In the dovetail ring of the optics carrier is a locating screw which must engage the notch in the accessory (page 33). This is best accomplished by moving the accessory slightly from side to side.
- ⇒ Release the clamping screw.
- ⇒ Introduce the accessory tube into the dovetail ring in the same manner as the binocular tube (page 33) and move it slightly from side to side until the locating screw engages the notch.
- ⇒ Tighten the clamping screw.
- ⇒ Fit the binocular tube on the accessory tube in the same manner.

### Phototubes

- ⇒ Fit the phototube as described above. The outlet can point to either the left or the right.
- Phototubes are designed to support modular photomicrographic systems, and video, film or SLR cameras.
  - The phototubes and the modular photomicrographic systems are detailed in a separate user manual, along with instructions for their assembly and use.
  - When you are not taking photographs, engage the observation beam path (VIS). When the photographic beam path (PHOT) is engaged, part of the light is directed to the camera, and the image seen down the binocular tube is darker.

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



## SLR camera

### Assembly

- 1 Phototube
- 2 Eyepiece tube
- 3 Photo eyepiece
- 4 Connecting tube
- 5 Screw the 0.32x camera objective, together with the appropriate T-mount adapter, to the adapter piece.
- 6 Secure the 35mm camera back on this assembly
- 7 Secure the combined unit (35mm camera back with assembly) on the 40mm connecting tube.

## TV- or film camera

### Assembly with photo eyepiece

- 1 Phototube
- 2 Eyepiece tube
- 3 10x or 16x photo eyepiece
- 4 Connect the f=50mm objective to the TV camera.
- 5 Secure the combined unit (TV camera with objective) on the phototube.

### Assembly without photo eyepiece

- 1 Phototube
- 2 Connect the TV camera to the C-mount adapter.
- 3 Insert the combined unit (TV camera with C-mount) into the phototube.

## LEICA 2, LEICA 3 video systems

- ⇒ Mount the video module on the optics carrier as described on page 48.

The individual system modules, their connections and their controls, are detailed in the separate user manual.

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## Accessories

### Drawing tube, discussion tube, filter-slide housing

#### Drawing tube

- The drawing tube is for both left-handers and right-handers. Fit the mirror on the left- or right-hand side of the optics carrier.
  - The drawing tube can only be used with the small incident-light stand and with the transmitted-light stand EB.
- ⇒ Fit the drawing tube between the optics carrier and the binocular tube as described on page 48.
- ⇒ Place a sheet of paper beneath the mirror.
- ⇒ Use a table lamp to illuminate the paper.

#### Drawing

- ⇒ Move the lever to the circle.
- ⇒ Look into the binocular tube, position the object, and focus on it.

- ⇒ Move the lever to the spot.
- ⇒ Look into the binocular tube and balance the illumination brightnesses so that the object, the drawing surface and the pencil are all seen together.
- In general, the paper needs to be illuminated strongly, but the object only weakly.
- ⇒ Move the lens mount of the drawing tube forwards or backwards until the drawing surface is sharply imaged.
- ⇒ Now you can draw the feature.
- If you place a graduated scale beneath the object and draw its intervals as well, you will easily determine the image scale on the paper.

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




### Discussion tube

Down the discussion tube two observers see the same erect, laterally-correct, three-dimensional image.

The assembly and use of the discussion tube are described in the separate user manual M2-261-004.

### Filter-slide housing

- The filter-slide housing can also be used beneath an accessory tube such as the phototube.
  - The filter-slide housing accepts strips of commercially-available gelatine filter, e.g. colour-compensation filters.
- ⇒ Secure the filter-slide housing as described on page 48.
- ⇒ Cut 23mmx45mm strips from the gelatine filters and insert them into the two filter holders of the filter slide.
- ⇒ Push the filter slide into its housing.

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

## Accessories

### Beam splitter

- Partition: 50% observation and 50% to the two lateral outlets
- The two lateral outlets can be used simultaneously and for different purposes.

### Assembly

- ⇒ Fit the beam splitter between the optics carrier and the binocular tube (page 48).
- ⇒ Release the clamping screw on the right or left outlet and remove the dust plug.
- ⇒ Carefully fit the accessory into the lateral outlet.
- The accessory is self-centring.

### Phototube

- ⇒ Secure the phototube to the left or right outlet of the beam splitter.
- ⇒ Attach the camera adapter to the SLR camera back (without camera objective).
- ⇒ Screw the f=250mm adapter (large fields) or the f=350mm camera adapter (high images scales) to the camera adapter.
- ⇒ Attach the camera and adapter together to the phototube and tighten the clamping screw.
- ⇒ Fit the format-indicating graticule for 35mm cameras into one of the eyepieces (page 54).

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



### Dual-station viewing

- ⇒ Attach the monocular tube for a second observer, or the stereo attachment with binocular tube, to the left or right outlet of the beam splitter.
- ⇒ Fit the eyepiece(s).
- ⇒ Adjust the dioptic setting(s) for the second observer.
- ⇒ Turn the second-observer tube as required.
- ⇒ Use the black rilled sleeve to match the image position to the direction of observation.

### Cine/TV tube

- ⇒ Attach the cine/TV tube to the left or right outlet of the beam splitter.
- ⇒ Screw the C-mount adapter to the cine or TV camera.
- ⇒ Attach the camera with adapter to the cine/TV tube and tighten the clamping screw.
- ⇒ Slacken the clamping screw to turn the outfit to any of four positions in accordance with requirements.
- ⇒ Retighten the clamping screw.

### Double-iris diaphragm

- The double-iris diaphragm increases the depth of field. It affects both the binocular tube and the accessories.
- ⇒ To increase the depth of field, move the lever from "OPEN" towards position 1.

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



## Miscellaneous accessories

### Attachment for vertical and oblique observation

#### Only for MS5, MZ6, MZ8

- The attachment for vertical and oblique observation enables spatial objects to be observed from above, and also from all round at an angle of 45°.
- The microscope carrier must be fitted in the lower position 1 (page 32).
- ⇒ For MS5, MZ6: Fit 1.0x achromatic objective by screwing its adapter ring directly into the microscope carrier.
- ⇒ To mount the attachment on the MZ8, two spacer rings are required.

The fitting and use of this accessory are described in detail in the user manual M2-265-004.

### Graticules

- The graticules for measuring and counting are in mounts which can be fitted directly into the eyepieces.
- A highly-precise stage micrometer with a 50mm scale and 0.1mm and 0.01 divisions enables the graticule to be calibrated in accordance with the magnification selected.
- It is best to measure with the microscope carrier AX in vertical position. The measurements are parallax-free and therefore more precise.
- Focusing and framing graticules indicate the format limits for all films.

The insertion of the graticules, the dioptic adjustment, and setting for parfocality when a graticule is in the eyepiece, are described in the user manuals for photomicrographic and measuring systems.

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## Stages

- Use the reducing ring diam. 120/80mm when fitting a stage into an incident-light stand, into the 20W transmitted-light stand, and into the transmitted-light stand for bright and dark field.

## Gliding stage

- The gliding stage can be moved in all directions and rotated completely.
- ⇒ Fit the gliding stage into the port in the base or into the reducing ring, and secure it with the eccentric lever.
- ⇒ Fit the 80mm diameter glass- or black/white stage plate, and attach the stage clips.

## Cup stage

- The cup stage can be inclined as required.
- You can secure objects to its rubberized upper surface with needles.
- The movable holder accepts petri dishes.
- ⇒ Fit the support ring into the port or into the reducing ring.
- ⇒ Place the cup stage on it.

## Polarization, transmitted light

- The analyzer, in combination with either the rotatable polarization stage or the polarizer on glass stage plate, enables double-refracting materials to be studied. These can include organic crystals, mineral grains, thin sections of rock or bone, polymers, stressed glass, and liquid crystals.
- Use the reducing ring diam. 120/80mm for the 20W transmitted-light stand and for the transmitted-light stand for bright and dark field.

The fitting and use of these accessories are described in detail in the user manual M2-265-004.

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
### Tips and hints What to do if ...

#### The field of view is partly shadowed.

##### Remedy

- Set the interpupillary distance correctly (page 16).
- Locate the eyepoint (page 17).

#### The image does not stay in focus

##### Remedy

- Fit the eyepieces correctly (page 33).
- Adjust your dioptric settings exactly as described (page 22).

#### The focusing drive drifts downwards or can only be moved with difficulty.

##### Remedy

- Adjust the ease of movement (pages 19, 37, 39).

#### If an electrical fault develops, first check the following:

- Is the voltage selector set correctly (bright-field transmitted-light stand, page 34; transformer, page 46)?
- Is the power switch on?
- Is the power cable attached correctly?
- Are all connecting cables attached correctly?
- Are the fuses in order?

#### The image is too dark.

##### Remedy

- Turn the regulating knob up high enough (bright-field transmitted light stand, page 24; transformer, page 46)?
- Open the double-iris diaphragm of the phototube.

### TV, photography

#### The image on the screen is too dark.

##### Remedy

- Adjust the video camera / monitor system in accordance with the directions of the manufacturer.
- Refer to the remedies listed above under: "The image is too dark".

### Unsharp photographs

##### Remedy

- Focus precisely (page 19).
- Adjust your dioptric settings exactly as described (page 22).

#### Colour photographs have a yellow cast.

##### Remedy

- Use artificial-light film.
- Turn the regulating knob of the transformer clockwise to the position of maximum voltage (bright-field transmitted-light stand, page 24; transformer, page 46).

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### Care

#### In this section

we will show you how to look after your microscope and how to clean it.

#### We guarantee the quality

You are working with a high-precision, high-performance instrument.

**We therefore guarantee the quality of each of our instruments. Our guarantee covers all faults in materials and manufacture. It does not, however, cover damage resulting from careless or improper handling.**

Please handle your valuable optical instrument with respect. It will then fulfil its function with the same precision for decades. Our instruments are famous for that!

If, contrary to expectations, your stereomicroscope develops a fault, please call in a specialist – contact your local Leica agent or the factory in Heerbrugg, Switzerland.

#### Protect your instruments

- from damp, vapours, acids, alkalis, and corrosive substances, and keep chemicals away from the instruments.
- from misuse. Never dismantle optical systems or mechanical parts unless instructions for doing so are given in the user manual.
- from oil and grease. Do not grease guide surfaces and mechanical parts.

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### **Dust and dirt will affect the quality of your results**

#### **Therefore**

- ⇒ Put a dust cover over the instrument during breaks in work.
- ⇒ Use dust caps to protect tube openings, tubes without eyepieces, and eyepieces.
- ⇒ Remove dust with a pneumatic rubber bulb or with a soft brush.
- ⇒ Clean objectives and eyepieces with special optic cleaning cloths and with pure alcohol.
- ⇒ When accessories are not being used, protect them against dust.

### **Cleaning polymer components**

Some components are made of polymer, or are polymer-coated. They are therefore pleasant and convenient to handle. The use of unsuitable cleaning agents and techniques can damage polymers, so always observe the following rules:

#### **Never clean polymers this way**

- Never clean in an ultrasonic bath. The polymer can become brittle and later break.
- Never use corrosive cleaning agents or ones containing acetone.
- Never use solvents for cleaning (except for ethanol and isopropanol).

#### **You can clean polymers harmlessly**

- with warm soapy water, wiping down afterwards with distilled water.
- with ethanol (industrial alcohol) and isopropanol.

When cleaning with ethanol or propanol, adhere to the appropriate safety regulations.

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**Notes**

**Total magnification and field diameter**

- V<sub>Ob</sub>** Magnification of objective
- V<sub>Oc</sub>** Magnification of eyepiece
- z** Magnification changer position
- t** Tube factor, e.g. 1.5x for coaxial incident light
- Sz** Field number of eyepiece.  
The field numbers are printed on the eyepieces:  
10x = 21, 16x = 14, 25x = 9.5, 40x = 6.

**Example**

- V<sub>Ob</sub>** 1.0x objective
- V<sub>Oc</sub>** 25x/9.5 eyepiece
- z** Zoom position 4.0
- t** Tube factor 1.5x (because of coaxial illuminator)

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
**Magnification down binocular tube:**

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